

Endothelin B receptor

Catalog No: #AP82855



Package Size: #AP82855-1 50ug #AP82855-2 100ug #AP82855-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelin B receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ET-B,Endothelin receptor non-selective type
Accession No.	Uniprot ID:P48302
Uniprot	P48302
GeneID	13618;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ednrb

Note: This product is for in vitro research use only