

Fatty acid-binding protein, intestinal

Catalog No: #AP82857



Package Size: #AP82857-1 50ug #AP82857-2 100ug #AP82857-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, intestinal
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Fatty acid-binding protein 2, Intestinal-type fatty acid-binding protein, I-FABP, FABPI
Accession No.	Uniprot ID:Q45KW7
Uniprot	Q45KW7
GeneID	595106;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FABP2

Note: This product is for in vitro research use only