

## Ectonucleoside triphosphate diphosphohydrolase 1

Catalog No: #AP82863



Package Size: #AP82863-1 50ug #AP82863-2 100ug #AP82863-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ectonucleoside triphosphate diphosphohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	NTPDase 1,Ecto-ATP diphosphohydrolase 1,Ecto-ATPDase 1,Ecto-apyrase,Lymphoid cell activation antigen,CD39,CD39
Accession No.	Uniprot ID:Q9MYU4
Uniprot	Q9MYU4
GeneID	397298;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ENTPD1

Note: This product is for in vitro research use only