

## DNA repair protein complementing XP-G cells

Catalog No: #AP82864



Package Size: #AP82864-1 50ug #AP82864-2 100ug #AP82864-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA repair protein complementing XP-G cells
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DNA excision repair protein ERCC-5,Xeroderma pigmentosum group G-complementing protein,ERCM2,XPG,XPGC
Accession No.	Uniprot ID:P28715
Uniprot	P28715
GeneID	2073;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ERCC5

Note: This product is for in vitro research use only