

Gastrotropin

Catalog No: #AP82865



Package Size: #AP82865-1 50ug #AP82865-2 100ug #AP82865-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastrotropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GT, 14 kDa bile acid-binding protein, Fatty acid-binding protein 6, I-BABP, Ileal lipid-binding protein, ILBP, Intestinal 15 kDa protein, I-15P, Ilbp, Illbp
Accession No.	Uniprot ID:P80020
Uniprot	P80020
GeneID	25440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fabp6

Note: This product is for in vitro research use only