

Endoplasmin

Catalog No: #AP82866



Package Size: #AP82866-1 50ug #AP82866-2 100ug #AP82866-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endoplasmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	94 kDa glucose-regulated protein,GRP-94,Endoplasmic reticulum resident protein 99,ERp99,Heat shock protein 90 kDa beta member 1,Polymorphic tumor rejection antigen 1,Tumor rejection antigen gp96,Grp94,Tra-1,Tra1
Accession No.	Uniprot ID:P08113
Uniprot	P08113
GeneID	22027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hsp90b1

Note: This product is for in vitro research use only