

Coagulation factor XII

Catalog No: #AP82867



Package Size: #AP82867-1 50ug #AP82867-2 100ug #AP82867-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Coagulation factor XII
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hageman factor, HAF
Accession No.	Uniprot ID:D3ZTE0
Uniprot	D3ZTE0
GeneID	306761;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F12

Note: This product is for in vitro research use only