

Endothelin-converting enzyme 1

Catalog No: #AP82878



Package Size: #AP82878-1 50ug #AP82878-2 100ug #AP82878-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelin-converting enzyme 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ECE-1
Accession No.	Uniprot ID:Q4PZA2
Uniprot	Q4PZA2
GeneID	230857;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ece1

Note: This product is for in vitro research use only