

Endothelin-1 receptor

Catalog No: #AP82879



Package Size: #AP82879-1 50ug #AP82879-2 100ug #AP82879-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endothelin-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Endothelin A receptor,ET-A,ET-AR
Accession No.	Uniprot ID:Q5KSU9
Uniprot	Q5KSU9
GeneID	450187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDNRA

Note: This product is for in vitro research use only