

Elongation factor 1-alpha 1

Catalog No: #AP82880

Package Size: #AP82880-1 50ug #AP82880-2 100ug #AP82880-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation factor 1-alpha 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	EF-1-alpha-1,Elongation factor Tu,EF-Tu,Eukaryotic elongation factor 1 A-1,eEF1A-1,EEF1A,EF1A
Accession No.	Uniprot ID:P68103
Uniprot	P68103
GeneID	282220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EEF1A1

Note: This product is for in vitro research use only