

Excitatory amino acid transporter 1

Catalog No: #AP82881

Package Size: #AP82881-1 50ug #AP82881-2 100ug #AP82881-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Excitatory amino acid transporter 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Glial high affinity glutamate transporter,High-affinity neuronal glutamate transporter,Sodium-dependent glutamate/aspartate transporter 1,GLAST-1,Solute carrier family 1 member 3,Eaat1,Gmt1
Accession No.	Uniprot ID:P56564
Uniprot	P56564
GeneID	20512;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc1a3

Note: This product is for in vitro research use only