

Fibroblast growth factor 4

Catalog No: #AP82884



Package Size: #AP82884-1 50ug #AP82884-2 100ug #AP82884-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	FGF-4,Heparin secretory-transforming protein,HST,Heparin-binding growth factor 4,HBGF-4,HST
Accession No.	Uniprot ID:P48803
Uniprot	P48803
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF4

Note: This product is for in vitro research use only