

Low affinity immunoglobulin gamma Fc region receptor II

Catalog No: #AP82886

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82886-1 50ug #AP82886-2 100ug #AP82886-3 1mg

Description

Product Name	Low affinity immunoglobulin gamma Fc region receptor II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fc gamma receptor IIB,Fc-gamma RII,IgG Fc receptor II beta,Lymphocyte antigen 17,Ly-17,CD32,Fcgr2b,Ly-17
Accession No.	Uniprot ID:P08101
Uniprot	P08101
GeneID	14130;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fcgr2

Note: This product is for in vitro research use only