

Dystonin

Catalog No: #AP82887



Package Size: #AP82887-1 50ug #AP82887-2 100ug #AP82887-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dystonin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	230 kDa bullous pemphigoid antigen,230/240 kDa bullous pemphigoid antigen,Bullous pemphigoid antigen 1,BPA,Dystonia musculorum protein,Hemidesmosomal plaque protein,BP230,BP240,BPAG1,DMH,DT,KIAA0728
Accession No.	Uniprot ID:Q03001
Uniprot	Q03001
GeneID	667;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DST

Note: This product is for in vitro research use only