

## Histone-lysine N-methyltransferase EZH2

Catalog No: #AP82888



Package Size: #AP82888-1 50ug #AP82888-2 100ug #AP82888-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histone-lysine N-methyltransferase EZH2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ENX-1,Enhancer of zeste homolog 2,Enx1h
Accession No.	Uniprot ID:Q61188
Uniprot	Q61188
GeneID	14056;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ezh2

Note: This product is for in vitro research use only