

N(G),N(G)-dimethylarginine dimethylaminohydrolase 1

Catalog No: #AP82890



Package Size: #AP82890-1 50ug #AP82890-2 100ug #AP82890-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	N(G),N(G)-dimethylarginine dimethylaminohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DDAH-1,DDAHI,Dimethylargininase-1,DDAH
Accession No.	Uniprot ID:O94760
Uniprot	O94760
GeneID	23576;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDAH1

Note: This product is for in vitro research use only