

Erythrocyte membrane protein band 4.2

Catalog No: #AP82892



Package Size: #AP82892-1 50ug #AP82892-2 100ug #AP82892-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Erythrocyte membrane protein band 4.2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Erythrocyte protein 4.2,E42P
Accession No.	Uniprot ID:P16452
Uniprot	P16452
GeneID	2038;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EPB42

Note: This product is for in vitro research use only