

## Endothelial protein C receptor

Catalog No: #AP82893



Package Size: #AP82893-1 50ug #AP82893-2 100ug #AP82893-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Endothelial protein C receptor   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Bovine   |
| Other Names        | Activated protein C receptor,APC receptor,Endothelial cell protein C receptor,CD201,EPCR |
| Accession No.      | Uniprot ID:Q28105  |
| Uniprot            | Q28105   |
| GeneID             | 282005;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PROCR

Note: This product is for in vitro research use only