

Separin

Catalog No: #AP82894



Package Size: #AP82894-1 50ug #AP82894-2 100ug #AP82894-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Separin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Caspase-like protein ESPL1,Extra spindle poles-like 1 protein,Separase,Esp1,Kiaa0165
Accession No.	Uniprot ID:P60330
Uniprot	P60330
GeneID	105988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Esp1

Note: This product is for in vitro research use only