

## Fatty acid-binding protein, adipocyte

Catalog No: #AP82895



Package Size: #AP82895-1 50ug #AP82895-2 100ug #AP82895-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fatty acid-binding protein, adipocyte
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Adipocyte lipid-binding protein,ALBP,Adipocyte-type fatty acid-binding protein,A-FABP,Fatty acid-binding protein 4
Accession No.	Uniprot ID:P48035
Uniprot	P48035
GeneID	281759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FABP4

Note: This product is for in vitro research use only