

## Endothelin-1 receptor

Catalog No: #AP82901



Package Size: #AP82901-1 50ug #AP82901-2 100ug #AP82901-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endothelin-1 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Endothelin A receptor,ET-A
Accession No.	Uniprot ID:Q29010
Uniprot	Q29010
GeneID	397457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EDNRA

Note: This product is for in vitro research use only