

Epsin-3

Catalog No: #AP82902



Package Size: #AP82902-1 50ug #AP82902-2 100ug #AP82902-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epsin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EPS-15-interacting protein 3
Accession No.	Uniprot ID:Q91W69
Uniprot	Q91W69
GeneID	71889;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Epn3

Note: This product is for in vitro research use only