

Protein delta homolog 1

Catalog No: #AP82903



Package Size: #AP82903-1 50ug #AP82903-2 100ug #AP82903-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein delta homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DLK-1,Adipocyte differentiation inhibitor protein,Preadipocyte factor 1,Pref-1,Dlk,Pref1,Scp-1
Accession No.	Uniprot ID:Q09163
Uniprot	Q09163
GeneID	13386;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dlk1

Note: This product is for in vitro research use only