

## Epithelial cell adhesion molecule

Catalog No: #AP82908



Package Size: #AP82908-1 50ug #AP82908-2 100ug #AP82908-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Epithelial cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Ep-CAM, Tumor-associated calcium signal transducer 1, TACSTD1, RCJMB04_29h4
Accession No.	Uniprot ID:Q5F381
Uniprot	Q5F381
GeneID	421292;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: EPCAM

Note: This product is for in vitro research use only