

Serine/threonine-protein kinase/endoribonuclease IRE1



Catalog No: #AP82910

Orders: order@signalwayantibody.com

Package Size: #AP82910-1 50ug #AP82910-2 100ug #AP82910-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase/endoribonuclease IRE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endoplasmic reticulum-to-nucleus signaling 1,Inositol-requiring protein 1,Ire1-alpha,IRE1a,Ire1
Accession No.	Uniprot ID:Q9EQY0
Uniprot	Q9EQY0
GeneID	78943;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ern1

Note: This product is for in vitro research use only