

Ferritin heavy chain

Catalog No: #AP82912



Package Size: #AP82912-1 50ug #AP82912-2 100ug #AP82912-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ferritin heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Ferritin H subunit,FTH
Accession No.	Uniprot ID:Q95MP7
Uniprot	Q95MP7
GeneID	100499480;403631;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FTH1

Note: This product is for in vitro research use only