

Glyceraldehyde-3-phosphate dehydrogenase

Catalog No: #AP82914



Package Size: #AP82914-1 50ug #AP82914-2 100ug #AP82914-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glyceraldehyde-3-phosphate dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GAPDH,Peptidyl-cysteine S-nitrosylase GAPDH,GAPD
Accession No.	Uniprot ID:P00356
Uniprot	P00356
GeneID	374193;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GAPDH

Note: This product is for in vitro research use only