

Ectonucleoside triphosphate diphosphohydrolase 1

Catalog No: #AP82916



Package Size: #AP82916-1 50ug #AP82916-2 100ug #AP82916-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ectonucleoside triphosphate diphosphohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	NTPDase 1,Ecto-ATP diphosphohydrolase 1,Ecto-ATPDase 1,Ecto-apyrase,Lymphoid cell activation antigen,CD39,CD39
Accession No.	Uniprot ID:O18956
Uniprot	O18956
GeneID	282223;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENTPD1

Note: This product is for in vitro research use only