

# Erythrocyte membrane protein band 4.2

Catalog No: #AP82917



Package Size: #AP82917-1 50ug #AP82917-2 100ug #AP82917-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Erythrocyte membrane protein band 4.2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Erythrocyte protein 4.2,BEP42
Accession No.	Uniprot ID:O46510
Uniprot	O46510
GeneID	281754;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:EPB42

Note: This product is for in vitro research use only