

FKHRL1 (Phospho-Ser644) Antibody

Catalog No: #12476

Package Size: #12476-1 50ul #12476-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FKHRL1 (Phospho-Ser644) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FKHRL1 (Phospho-Ser644) Antibody detects endogenous levels of FKHRL1 only when phosphorylated at Ser644
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FKHRL1 (Phospho-Ser644)
Target Name	FKHRL1
Modification	Phospho
Other Names	FOXO3, AF6q21, AF6q21 protein, Forkhead box O3, Forkhead box O3A, FOXO3A, FKHRL1P2, Forkhead box protein O3, FOXO2, FKHRL1
Accession No.	Swiss-Prot#: O43524NCBI Gene ID: 2309
Uniprot	O43524
GeneID	2309;
Target Species	human
Calculated MW	97kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of 1. Hela cell and 2. Hela cell treated with 100ng/mL LPS for 30min, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only