

Endoplasmic reticulum metallopeptidase 1

Catalog No: #AP82921



Package Size: #AP82921-1 50ug #AP82921-2 100ug #AP82921-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endoplasmic reticulum metallopeptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Felix-ina,FXNA,KIAA1815
Accession No.	Uniprot ID:Q7Z2K6
Uniprot	Q7Z2K6
GeneID	79956;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ERMP1

Note: This product is for in vitro research use only