

## Receptor tyrosine-protein kinase erbB-4

Catalog No: #AP82924



Package Size: #AP82924-1 50ug #AP82924-2 100ug #AP82924-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Receptor tyrosine-protein kinase erbB-4  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Proto-oncogene-like protein c-ErbB-4, Tyrosine kinase-type cell surface receptor HER4, p180erbB4, HER4 |
| Accession No.      | Uniprot ID:Q15303  |
| Uniprot            | Q15303   |
| GeneID             | 2066;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ERBB4

Note: This product is for in vitro research use only