

# Adropin

Catalog No: #AP82925



Package Size: #AP82925-1 50ug #AP82925-2 100ug #AP82925-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Adropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Energy homeostasis-associated protein
Accession No.	Uniprot ID:A2VE22
Uniprot	A2VE22
GeneID	783487;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: ENHO

Note: This product is for in vitro research use only