

Low affinity immunoglobulin gamma Fc region receptor III

Catalog No: #AP82926

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP82926-1 50ug #AP82926-2 100ug #AP82926-3 1mg

Description

Product Name	Low affinity immunoglobulin gamma Fc region receptor III
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IgG Fc receptor III,Cytolytic trigger molecule G7,Fc-gamma RIII,FcRIII,CD16
Accession No.	Uniprot ID:Q28942
Uniprot	Q28942
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCGR3

Note: This product is for in vitro research use only