

Proto-oncogene c-Fos

Catalog No: #AP82932



Package Size: #AP82932-1 50ug #AP82932-2 100ug #AP82932-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene c-Fos
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Cellular oncogene fos
Accession No.	Uniprot ID:O97930
Uniprot	O97930
GeneID	100144486;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOS

Note: This product is for in vitro research use only