

Ephrin type-A receptor 3

Catalog No: #AP82939



Package Size: #AP82939-1 50ug #AP82939-2 100ug #AP82939-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ephrin type-A receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EPH-like kinase 4,EK4,Tyrosine-protein kinase TYRO4,Tyrosine-protein kinase receptor ETK1,Etk1,Mek4,Tyro4
Accession No.	Uniprot ID:P29319
Uniprot	P29319
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Epha3

Note: This product is for in vitro research use only