

Beta-defensin 1

Catalog No: #AP82941



Package Size: #AP82941-1 50ug #AP82941-2 100ug #AP82941-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	BD-1,Defensin, beta 1
Accession No.	Uniprot ID:O89117
Uniprot	O89117
GeneID	83687;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Defb1

Note: This product is for in vitro research use only