

Forkhead box protein C2

Catalog No: #AP82945

Package Size: #AP82945-1 50ug #AP82945-2 100ug #AP82945-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein C2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Forkhead-related protein FKHL14,Mesenchyme fork head protein 1,MFH-1 protein,Transcription factor FKH-14,FKHL14,MFH1
Accession No.	Uniprot ID:Q99958
Uniprot	Q99958
GeneID	2303;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOXC2

Note: This product is for in vitro research use only