

Formyl peptide receptor 2

Catalog No: #AP82946

Package Size: #AP82946-1 50ug #AP82946-2 100ug #AP82946-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Formyl peptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Formylpeptide receptor-related sequence 2,Lipoxin A4 receptor-like protein,N-formylpeptide receptor-like 2,Fpr-rs2
Accession No.	Uniprot ID:O88536
Uniprot	O88536
GeneID	14289;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fpr2

Note: This product is for in vitro research use only