

Fibrinogen-like protein 1

Catalog No: #AP82948

Package Size: #AP82948-1 50ug #AP82948-2 100ug #AP82948-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibrinogen-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HP-041,Hepassocin,Hepatocyte-derived fibrinogen-related protein 1,HFREP-1,Liver fibrinogen-related protein 1,LFIRE-1,HFREP1
Accession No.	Uniprot ID:Q08830
Uniprot	Q08830
GeneID	2267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGL1

Note: This product is for in vitro research use only