

Focal adhesion kinase 1

Catalog No: #AP82949

Package Size: #AP82949-1 50ug #AP82949-2 100ug #AP82949-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Focal adhesion kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FADK 1,Focal adhesion kinase-related nonkinase,FRNK,Protein-tyrosine kinase 2,p125FAK,pp125FAK,Fadk,Fak,Fak1,Kiaa4203
Accession No.	Uniprot ID:P34152
Uniprot	P34152
GeneID	14083;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptk2

Note: This product is for in vitro research use only