

Coagulation factor XIII A chain

Catalog No: #AP82951



Package Size: #AP82951-1 50ug #AP82951-2 100ug #AP82951-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coagulation factor XIII A chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Coagulation factor XIIIa, Protein-glutamine gamma-glutamyltransferase A chain, Transglutaminase A chain, F13A
Accession No.	Uniprot ID:P00488
Uniprot	P00488
GeneID	2162;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: F13A1

Note: This product is for in vitro research use only