

Histone-lysine N-methyltransferase EZH2

Catalog No: #AP82953



Package Size: #AP82953-1 50ug #AP82953-2 100ug #AP82953-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone-lysine N-methyltransferase EZH2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ENX-1,Enhancer of zeste homolog 2,Lysine N-methyltransferase 6,KMT6
Accession No.	Uniprot ID:Q15910
Uniprot	Q15910
GeneID	2146;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EZH2

Note: This product is for in vitro research use only