

Low affinity immunoglobulin gamma Fc region receptor II

Catalog No: #AP82955

Package Size: #AP82955-1 50ug #AP82955-2 100ug #AP82955-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low affinity immunoglobulin gamma Fc region receptor II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IgG Fc receptor II,Fc-gamma RII,FcRII,CD32
Accession No.	Uniprot ID:Q28110
Uniprot	Q28110
GeneID	282229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCGR2

Note: This product is for in vitro research use only