

Fibroblast growth factor 5

Catalog No: #AP82956



Package Size: #AP82956-1 50ug #AP82956-2 100ug #AP82956-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-5,Heparin-binding growth factor 5,HBGF-5,Smag-82
Accession No.	Uniprot ID:P12034
Uniprot	P12034
GeneID	2250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF5

Note: This product is for in vitro research use only