

Fibroblast growth factor-binding protein 1

Catalog No: #AP82962



Package Size: #AP82962-1 50ug #AP82962-2 100ug #AP82962-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGF-BP,17 kDa heparin-binding growth factor-binding protein,17 kDa HBGF-binding protein,FGFBP,HBP17
Accession No.	Uniprot ID:Q14512
Uniprot	Q14512
GeneID	9982;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGFBP1

Note: This product is for in vitro research use only