

Low affinity immunoglobulin epsilon Fc receptor

Catalog No: #AP82964



Package Size: #AP82964-1 50ug #AP82964-2 100ug #AP82964-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low affinity immunoglobulin epsilon Fc receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fc-epsilon-RII,Lymphocyte IgE receptor,CD23,Fcεr2a
Accession No.	Uniprot ID:P20693
Uniprot	P20693
GeneID	14128;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fcεr2

Note: This product is for in vitro research use only