

Dual specificity protein phosphatase 1

Catalog No: #AP82965



Package Size: #AP82965-1 50ug #AP82965-2 100ug #AP82965-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity protein phosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mitogen-activated protein kinase phosphatase 1,MAP kinase phosphatase 1,Protein-tyrosine phosphatase CL100,Protein-tyrosine phosphatase non-receptor type 16,Cl100,Mkp1,Ptpn16
Accession No.	Uniprot ID:Q64623
Uniprot	Q64623
GeneID	114856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dusp1

Note: This product is for in vitro research use only