

# Ectonucleotide pyrophosphatase/phosphodiesterase family member 2

Catalog No: #AP82966

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP82966-1 50ug #AP82966-2 100ug #AP82966-3 1mg

## Description

Product Name	Ectonucleotide pyrophosphatase/phosphodiesterase family member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	E-NPP 2,Autotaxin,Extracellular lysophospholipase D,LysoPLD
Accession No.	Uniprot ID:A1A4K5
Uniprot	A1A4K5
GeneID	532663;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ENPP2

Note: This product is for in vitro research use only