

Ectonucleoside triphosphate diphosphohydrolase 5

Catalog No: #AP82968



Package Size: #AP82968-1 50ug #AP82968-2 100ug #AP82968-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ectonucleoside triphosphate diphosphohydrolase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NTPDase 5,CD39 antigen-like 4,ER-UDPase,Guanosine-diphosphatase ENTPD5,GDPase ENTPD5,Nucleoside diphosphatase,Uridine-diphosphatase ENTPD5,UDPase ENTPD5,CD39L4,PCPH
Accession No.	Uniprot ID:O75356
Uniprot	O75356
GeneID	957;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENTPD5

Note: This product is for in vitro research use only